


MANDREL ORDERING INSTRUCTIONS



Our plugs come with the standard mandrels listed below. To order a different mandrel configuration, please follow these instructions:

add “-SM” to Item for small mandrel

add “-OT” to Item for stubby mandrel

add “ -1” to Item to add a 1” thru hole to our standard mandrel

New orders receive free mandrel retrofit

BK PLUGS - standard mandrel

ITEM	HEIGHT	TOP DIAMETER	BOTTOM DIAMETER
BK1-12	12"	5.1"	8"
BK2-14	14"	6"	10.25"
BK1	15"	4.5"	8"
BK2	18"	4.5"	10.25"

SPINAL PLUGS - mandrel with 1" thru hole

ITEM	HEIGHT	AP	ML
		TOP - BOTTOM	TOP - BOTTOM
T1	20"	9"	10"
T2	23"	12"	13"
T2T	28"	12"	13"
T3	28"	14"	15"
T3H*	28"	16.5"	17.5"
T4	30.5"	22.5	22.5

*T3H - Choose for Pendulous Abdomen

KNEE/SYMES PLUGS - standard mandrel

ITEM	HEIGHT	TOP DIAMETER	BOTTOM DIAMETER
K1	20"	6"	8"
K2	20"	7.5"	10"
K3	22.5"	8.5"	13"
SYME1	24"	6"	8.25"

CRANIAL PLUGS - small mandrel

ITEM	HEIGHT	TOP DIAMETER	BOTTOM DIAMETER
CT8	8"	8"	8"
CT10	10"	10"	10"

AK PLUGS - standard mandrel

ITEM	HEIGHT	TOP DIAMETER	BOTTOM DIAMETER
CM1	16.5"	7"	11.5"
CM2-18	18"	8"	13"
CM2	22.5"	8"	13"
CM2-E	24"	8"	13"

AFO PLUGS - standard mandrel

ITEM	HEIGHT	TOP DIAMETER	BOTTOM DIAMETER
AFO1	30"	6.5"	8.25"
AFO2	36"	7"	9"

CUSTOM SIZES & CONFIGURATIONS AVAILABLE

BioSculptor®
innovative solutions